



Rolence Enterprise, Inc.  
% Sterling Cheng  
General Manager  
18-3 Lane 231 Pu Chung Rd., Chungli  
Taoyuan, 32083  
TAIWAN

March 16, 2018

Re: K172928  
Trade/Device Name: XR-01 Portable X-ray System  
Regulation Number: 21 CFR 872.1800  
Regulation Name: Extraoral source x-ray system  
Regulatory Class: II  
Product Code: EHD  
Dated: February 22, 2018  
Received: February 23, 2018

Dear Sterling Cheng:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Robert Ochs, Ph.D.

Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K172928

Device Name

XR-01 Portable X-Ray System

Indications for Use (Describe)

The XR-01 Portable X-Ray System is intended to be used by trained dentists and dental technicians as an extraoral x-ray source for producing diagnostic x-ray images using intraoral image receptors. Its use is intended for both adult and child subjects

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

**Date 510(k) Summary prepared** 2018/2/1  
**Submitter Name** Rolence Enterprise Inc.  
 Address: 18-3 Lane 231 Pu Chung Rd., Chungli, Taoyuan,  
 Taiwan  
**Contact Person** Jui-Ya Lai  
 Regulatory Affairs Quality Assurance Associate Manager  
**Telephone** +886-3-4631999  
**FAX** +886-3-4631997

### Device Information

**Type of 510(k) Submission** Traditional  
**Trade or Proprietary Name** XR-01 Portable X-Ray System  
**Common or Usual Name** Portable X-ray System  
**Regulation Name** Extraoral source x-ray system  
**Regulation Number** 872.1800  
**Product Code** EHD  
**Class of Device** Class II  
**Panel** Radiology

### 1. Comparison with predicate devices:

	<b>Subject Device</b>	<b>Primary Predicate Device</b>	<b>Reference Devices</b>		
<b>Device Name</b>	XR-01 Portable X-Ray System	Nomad Pro Intraoral X-Ray system	Port-X II	DX-3000 Portable X-Ray System	
<b>Regulation Number</b>	872.1800	872.1800	872.1800	872.1800	
<b>Regulation Name</b>	Extraoral source x-ray system	Extraoral source x-ray system	Extraoral source x-ray system	Extraoral source x-ray system	
<b>Regulatory Class</b>	II	II	II	II	
<b>Primary Product Code</b>	EHD	EHD	EHD	MUH	
<b>Applicant Name</b>	Rolence Enterprise Inc.	Aribex Inc.	Genoray Co., Ltd.	Dexcowin Co., Ltd.	
<b>510(k) Number</b>	-	K081664	K063121	K133007	
<b>INTENDED USE</b>	All four systems are intended as extraoral X-ray sources to be used with intraoral image receptors for diagnostic imaging by dentist or dental technicians.				
<b>Mechanical</b>	<b>Size (mm)</b>	174 x 178 x 257	266.7 x 133.35 x 304.8	197 x 147 x 145.5	139 x 163.5 x 66.5

	<b>Source to skin distance</b>	200 mm	210 mm	100 mm	> 100 mm (fixed by cone)
	<b>X-ray field Size</b>	60 mm round	60 mm round	70 mm round	50 mm round
	<b>User Interface</b>	Up-down buttons for exposure time selection with timer display. Additionally, several user-selectable preset times with patient size, image receptor type and tooth selection icon on control panel.	Up-down buttons for exposure time selection, with timer display. Additionally, several user-selectable preset times with patient size, image receptor type, and tooth selection icons on an LCD display.	Up-down buttons for exposure time selection with timer display. Additionally, several user-selectable preset times with patient size and tooth selection icons on an LCD display.	Up-down buttons for exposure time selection with timer display. Additionally, several user-selectable preset times with patient size and tooth selection icons on an LCD display.
	<b>Backscatter radiation protection</b>	Available as optional accessory	171.45 mm dia. Pb-filled acrylic plastic scatter shield	None	Available as optional accessory
	<b>Exposure Switch</b>	Exposure button at front body	Trigger on Handset	Exposure button at x-ray control panel	Exposure button at front body
	<b>Tube head mounting</b>	Handheld or on a tripod	Handheld	Handheld or on a tripod	Handheld
<b>Electrical</b>	<b>Energy source</b>	Rechargeable 25.2 V DC Lithium Polymer battery pack	Rechargeable 22.2 V DC Lithium Polymer battery pack	Rechargeable 22.2 V DC Lithium Polymer battery pack	Rechargeable 16.8 V DC Lithium Polymer battery pack
	<b>Exposure time</b>	0.01 – 2.00 seconds in 0.01 steps	0.02 - 1.00 seconds in 0.01 increments	0.01 – 2.0 seconds in 46 steps	0.05 – 1.35 seconds in 0.01 steps
	<b>mA</b>	2 mA fixed	2.5 mA fixed	2 mA fixed	2 mA
	<b>kVp</b>	60 kVp fixed	60 kVp fixed	60 kVp fixed	60 kvp
	<b>Waveform</b>	Constant Potential (DC)	Constant Potential (DC)	Constant Potential (DC)	Constant Potential (DC)
	<b>Duty Cycle</b>	1:30	1:60	N/A	1:60
<b>Applied Standards</b>		IEC 60601-1,	IEC/EN	UL/IEC	21 CFR

	EN 60601-1-2, IEC 60601-1-3, IEC 60601-2-65, IEC 62133 21 CFR 1020.30, 21 CFR 1020.31	60601-1, EN 60601-1-2, IEC 60601-1-3, IEC 60601-2-7, 21 CFR 1020.30, 21 CFR 1020.31	60601-1, IEC 60601-1-2 IEC 60601-2-7, 21 CFR 1020.30, 21 CFR 1020.31	1020.30, 21 CFR 1020.31
--	--	--	--	-------------------------------

### Device Description

XR-01, a portable X-ray system is mainly designed for dental diagnosis via X-ray exposure with intraoral image receptors, supplied with 25.2 V DC rechargeable Li-ion polymer battery. The portable X-ray system consist of an x-ray generator with an x-ray tube, device controller, beam limiting device, AC-DC power supply, charging unit and optional back scatter shield.

### Indication for Use

The XR-01 Portable X-Ray System is intended to be used by trained dentists and dental technicians as an extraoral x-ray source for producing diagnostic x-ray images using intraoral image receptors. Its use is intended for both adult and child subjects.

### Performance Testing

Subject device XR-01 complies with safety and performance standards which are in accordance with 21 CFR 1020.30 and 31. Performance test validate that accuracy, reproducibility, deviation and linearity of technical factors such as tube voltage, tube current and exposure time is meet the criteria. Safety tests in terms of HVL, leakage, x-ray field size meet Federal standard requirements.

### Standards and Guidance Documents

#### Electrical Safety and EMC Standards

- IEC 60601-1:2005 (Third Edition) + CORR. 1:2006 + CORR. 2:2007 + A1:2012
- EN 60601-1-2:2015
- IEC 60601-1-3:2008 (Second Edition) + A1:2013 for use in conjunction with IEC 60601-1:2005 (Third Edition) + A1:2012
- IEC 60601-1-6:2010 (Third Edition) + A1:2013
- IEC 60601-2-65:2012 (First Edition) for use in conjunction with IEC 60601-1:2005 (Third edition)

#### Battery Safety Standards

- IEC 62133:2012 (Second Edition)

#### Usability Standards

- IEC 62366: 2007 (First Edition) + A1: 2014

#### Radiation Control

- 21 CFR 1020.30
- 21 CFR 1020.31

**Guidance**

- Radiation Safety Considerations for X-Ray Equipment Designed for Hand-Held Use
- How to Prepare a Special 510(k)
- Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices
- Pediatric Information for X-ray Imaging Device Premarket Notifications
- Applying Human Factors and Usability Engineering to Medical Devices
- Content of Premarket Submissions for Management of Cybersecurity in Medical Devices
- Refuse to Accept Policy for 510(k)s
- eCopy Program for Medical Device Submissions

**Conclusion**

The subject device and predicate devices have same intended use and principle of operation. Also, both of them met Federal regulations 21 CFR Part 1020.30, 1020.31 and international medical device safety standard. In accordance with the Federal Food, Drug and Cosmetic Act, 21 CFR Part 807 and based on the information provided in this premarket notification Rolence Enterprise Inc. concludes that the Portable X-Ray System (XR-01) is safe and effective and substantially equivalent to predicate devices as described herein.